

**DECK MOUNT**

Saving millions of lives



 **respirare**  
sleep magic

**Respirare VT 70D** is designed for people who want to enjoy high-quality sleep any-time, anywhere. It's compact and light weight design, making it easy to carry.



# Respirare VT 70D

Sleep Better, Anywhere, Everywhere!



**Respirare VT 70D** is a compact and light weight device, making it easily portable and convenient to carry.

Respirare VT 70D is equipped with advanced sensors that can detect even the slightest changes in breathing pattern. It has been rigorously tested and calibrated to ensure accurate and reliable monitoring of breathing events including snoring.

### Humidification (HEM Technology)

Humidification without Electricity and Maintenance . A cutting edge humidification system, it enables the device to provide effective humidification without requiring water or electricity.

## Product Code

Respirare Auto CPAP

PAP VT 70 D

Note: Technical specification may change without notice. Always refer to your device user manual.

## Technical Specification

Modes	CPAP, Auto CPAP
Operating Pressure Range	4-20 cm H2O
Tubing Length	6 ft 5 in
Tubing Inner Diameter	15 mm
Sound - Pressure Level	30 dBA
Sound Level	37 dBA
Machine Dimensions	14.3 cm x 9.0 cm x 5.4 cm
Machine Weight	317 gms
Power Supply Weight	45 lbs (7.2 oz)
AC Input Range	220V 50 Hz
Typical Electricity Consumption	10 W
Onboard Storage Memory	1 Year
Standby Electricity Consumption	0.3 W
Operating Temperature	+41°F to +95°F (+5°C to +35°C)
Operating Humidity	10 to 93% relative Humidity, non- condensing
Operating Altitude	Sea level to 8,500 ft
Storage and Transport Temp.	-13°F to +158°F (-25°C to +70°C)
Storage and Transport Humidity	5 to 95% relative humidity, non condensing
Repairable	Yes
Warranty	2 Years
Conforms to Standards	IEC 60601-1-1 IEC 60601-1-2 IEC 80601-2-12
	Class 2 (double insulated, Type BF, Ingress Protection IP22)

\*Term & Conditions Apply

- Ensures a peaceful sleep environment
- Real Time Data Tracking and Analysis.
- Similar to i-mode/ A-Flex
- iOS and Android software operated



53.3g  
Lightweight



Secure Fit



Flexibility

